

New England Biolabs Product Specification

Product Name: Oligo d(T)₂₅ Magnetic Beads
Catalog #: S1419S
Concentration: 5 mg/ml
Shelf Life: 36 months
Storage Temp: 4°C
Storage Conditions: 0.02 % NaN₃, 0.05 % Tween®20, 1 X PBS, (pH 7.4 @ 25°C)
Specification Version: PS-S1419S v2.0
Effective Date: 25 Oct 2021

Assay Name/Specification (minimum release criteria)

Binding Capacity (Magnetic Beads) - Oligo d(T)₂₅ Magnetic Beads (500 µg) were equilibrated and incubated with 100 µl of 67 µg/ml rA30 for 15 minutes at 25°C, then washed and the rA30 eluted. Binding capacity was determined to be >5 µg of rA30 per mg of beads.

Functional Testing (mRNA Isolation) - Oligo d(T)₂₅ Magnetic Beads were equilibrated and incubated with freshly prepared eukaryotic cell lysate for direct mRNA isolation. The beads were washed and the mRNA eluted. The eluate was evaluated on an Agilent Bioanalyzer and the enriched poly(A)+ RNA contains ≤15% rRNA.

Non-Specific DNase Activity (16 hour, Buffer) - A 50 µl reaction in Oligo d(T)₂₅ Magnetic Bead Storage Buffer containing 1 µg of PhiX174-HaeIII DNA incubated for 16 hours at 37°C results in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis.

RNase Activity (Extended Digestion) - A 10 µl reaction in NEBuffer 4 containing 40 ng of a 300 base single-stranded RNA and a minimum of 1 µl of Oligo d(T)₂₅ Magnetic Beads is incubated at 37°C. After incubation for 16 hours, >90% of the substrate RNA remains intact as determined by gel electrophoresis using fluorescent detection.

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 25 Oct 2021

Derek Robinson
Director, Quality Control

